

Description	Measurements ca. (mm)	Water volume ca. (ccm)	Weight (only cuff) ca. (kg)	Weight (incl. Accessories) ca. (kg)	Part-N° (incl. Accessories)	EAN-Code (incl. Accessories)
Eye cuff - single use - incl. 1 Hilofix 25	230 x 90	29	0,018	0,028	15000002	42 503871 0046 6
Eye cuff, open - single use - incl. 1 Hilofix 25	225 x 95	41	0,018	0,028	15000003	42 503871 0047 3
Small round cuff - single use - incl. 2 Hilofix 25	Ø 200	49	0,029	0,046	15000007	42 503871 0048 0
Round cuff - single use - incl. 2 Hilofix 25	Ø 240	90	0,039	0,056	15000008	42 503871 0049 7
Large surface cuff -single use - incl. 2 Hilofix 75	700 x 250	377	0,105	0,159	15000009	42 503871 0052 7
Small surface cuff - single use - incl. 1 Hilofix 25	240 x 70	33	0,031	0,041	15000010	42 503871 0053 4
Lower face cuff with link - single use - incl. 2 Hilofix 25	350 x 220	74	0,031	0,048	15000012	42 503871 0054 1
Upper face cuff - single use - incl. 3 Hilofix 25	340 x 150	60	0,029	0,052	15000013	42 503871 0055 8
Lower leg cuff - single use - incl. 2 Hilofix 75	500 x 340	349	0,104	0,158	15000015	42 503871 0056 5
Knee cuff - single use - incl. 2 Hilofix 75	470 x 300	241	0,076	0,130	15000016	42 503871 0057 2
Foot cuff - single use - incl. 2 Hilofix 75	470 x 360	331	0,084	0,138	15000018	42 503871 0058 9
Nose cuff - single use - incl. 2 Hilofix 25	235 x 120	20	0,021	0,038	15000019	42 503871 0059 6
Nose cuff T-cast - single use - incl. 1 Hilofix 25	270 x 160	66	0,023	0,033	15000020	42 503871 0060 2
Lower face cuff without link - single use - incl. 2 Hilofix 25	310 x 190	56	0,030	0,047	15000022	42 503871 0061 9
Universal cuff mini -single use - incl. 1 Hilofix 25/65	100 x 100	17	0,013	0,023	15000004	42 203871 0113 2
Trapezoidal cuff - single use - incl. 2 Hilofix 25	420 X 320	112	0,083	0,135	15000014	43 503871 0117 0